

TravelPlus

20-MINUTES SESSION



**WHO DOESN'T WANT TO MISS OUT THIS ENERGY-THERAPY
WHEN BEING ON A JOURNEY CAN USE THE TRAVEL PLUS
ANYWHERE AT ANY TIME TO GET NEW ENERGY.**

The Airnergy Travel Plus is equipped with 4 AIRNERGY activation units and colour light therapy (blue) and is battery-powered.
The recommended application time is 20 minutes.

EQUIPMENT

4 AIRNERGY catalysts (activation units)
Integrated colour light therapy (blue)
3 intensity levels (50 %, 75 %, 100 %)
LCD display
Five-language menu guidance (DE, EN, FR, IT, ES)
Microcontroller-controlled electronics
Individual programme settings
Automatic self-test
Switchable acoustic signal
LED strip (battery indicator)
Main switch
Battery powered
Optional: aromatherapy

TECHNICAL SPECIFICATIONS

Power supply	100-240 V
Device battery	built-in lithium-polymer battery 3s/5,800
Battery power	approx. 210 min.
Power consumption	11.7 Volt, max. 5,800 mA
Air supply	approx. 4 l per min.
Dimensions (mm) case	386 x 245 x 240 (W x H x D)
Dimensions (mm) battery charger	80 x 160 x 60 (W x H x D)
Battery charger	charge end voltage: 12.6 V
Weight	approx. 5.25 kg

INFORMATION

SERVICE	Every 2,000 hours of operation (approx. 6,000 sessions)
GUARANTEE	2 years (no guarantee on battery)

SHIPMENT INCLUDES 1 battery charger, 20 nasal cannulas, 1 glass bottle incl. sparkling element, 2 air filters, 1 cannula bag, 1 user manual

ACCESSORIES (OPTIONAL) Aroma bottle

ARTICLE-NO. 2090055



AIRNERGY IRELAND

AIRNERGY+
ENERGY FROM INSIDE

PH:01-8206346 FAX:01-8206347 EMAIL: info@airnergyireland.ie WEBSITE:www.airnergyireland.ie